Airnet® IIs

Stainless Steel 2 Channel Particle Sensor



The Airnet® II_s particle sensor makes it easy and cost-effective to monitor your cleanroom. This particle sensor offers a small footprint, unparalleled performance, and data transmission capabilities while meeting the specification of ISO 21501-4 and ISO 14644-1:2015, all enclosed in a 316L Stainless Steel case.

Simple installation with one single cable, using Power Over Ethernet (PoE) Technology. Communication capabilities include Ethernet to interface with Pharmaceutical Net, Facility Net, or FacilityPro® software, OPC communications, and Modbus TCP/IP communications.

Data integrity is maintained through the use of a data queue feature that continues to gather data even if network communication is lost.

To ensure proper flow conditions and vacuum system operation, these units incorporate a Dynamic Flow Sensing system that will alarm with a 15% change in flow conditions.

For applications where decontamination using Vaporized Hydrogen Peroxide (VHP) is required, an optional VHP-compatible unit is available for a simple installation without complex valving (Airnet 510XR).

BENEFITS

- Proven technology provides reliable and accurate data
- Allows for immediate reaction to particle contamination events
- A low-cost solution for multipoint monitoring
- Interfaces with Facility Net, Pharmaceutical Net and FacilityPro Software for comprehensive management of cleanroom conditions
- A small footprint and wall mounting bracket make it easy to install in cleanrooms and mini-environments.
- A laser diode (LD) drastically reduces the need for maintenance and extends product lifecycle.
- Automatic laser shutdown reduces laser failures.
- Data queue maintains data integrity when communication is lost.
- Optional XR coating protects sensors against corrosive or oxidizing vapors in VHP sterilization processes.

FEATURES

- 2 channels
- 0.5 and 5.0 µm size range, suitable for Pharmaceutical application
- 1.0 CFM flow rate
- Interfaces with object linking and embedding (OLE) for process control (OPC) communications, Modbus communications
- Chemical-resistant 316L Stainless Steel enclosure
- Low sample point cost
- Small enough for use in remote locations
- ISO 14644-1:2015 Compliant
- ISO 21501-4 Compliant

APPLICATIONS

- Cleanroom monitoring
- Dedicated monitoring of critical locations
- Trend analysis
- Statistical process control
- · Multi-location monitoring
- Isolator monitoring



Without measurement there is no control

Airnet® II_s

Stainless Steel 2 Channel Particle Sensor

1

| | 510s | 510s XR |
|---------------------------------------|---|------------------------------------|
| Size Range | 0.5, 5.0 μm | |
| Flow Rate | 1.0 CFM (28.3 LPM) | |
| Counting Efficiency | 50% ± 20% for most-sensitive channel. Meets ISO 21501-4 $100%$ ± 10% at 1.5 to 2.0 times channel one size. Meets ISO 21501-4. | |
| Zero Count | ≤ 7.07 counts/m³ | |
| Maximum Concentration ¹ | 957,824/ft³ | |
| Laser Source | Diode | |
| Laser Classification | Class 1 per EN60825 (Internally, a Class IIIB laser is used, per EN60825) | |
| Exterior Surface | 316L Stainless Steel | |
| Dimensions (H x W x L) | 3.5 x 7.1 x 4.0 in (8.9 x 18.0 x 10.2 cm) without ISP and bracket | |
| Weight | 2.4 lb (1.08 kg) | |
| Vacuum Connection | 8 mm OD tubing | |
| Flow System | External vacuum required; automatic laser shutoff and alarm on 15% flow variation | |
| Vacuum Source | > 15 in Hg (> 410 mBar) below atm pressure | |
| Power | Power over Ethernet 48 VDC via a PoE router | |
| Communication Connectors | Modbus TCP) | Systems proprietary protocol, OPC, |
| Status Indicators | RS-232 (configuration and d | |
| Calibration | Programmable status (two-color LED), Activity (one-color LED) Calibration materials used are traceable to the National Institute | |
| Calibration | | gy (NIST) and meet ISO 21501-4 |
| Environment | Temperature: 4 – 35 °C (32 – 95 °F) Humidity: 5 – 95% non-condensing relative | |
| Complies with | EU RoHS, ISO 21501-4, ISO 14644-1:2015 | |



Airnet II with wall mounting bracket



316L Stainless Steel Wall Mounting Bracket, included

¹ 10% coincidence loss at maximum concentration. Airnet® and FacilityPro® are registered trademarks of Particle Measuring Systems, Inc. All other trademarks are the property of their respective owners Particle Measuring Systems, Inc. reserves the right to change specifications without notice. © 2016 Particle Measuring Systems, Inc.

HEADQUARTERS

5475 Airport Blvd Boulder, Colorado 80301 USA T: +1 303 443 7100, +1 800 238 1801

Instrument Service & Support T: +1 800 557 6363

Customer Response Center T: +1 877 475 3317

E: info@pmeasuring.com

www.pmeasuring.com info@pmeasuring.com

GLOBAL OFFICES

AUSTRIA

T: +43 512 390 500

E: pmsaustria@pmeasuring.com

BENELUX

T: +32 10 23 71 56

E: pmsbelgium@pmeasuring.com

T: +55 11 5188 8227

E: pmsbrazil@pmeasuring.com

CHINA

T: +86 21 6113 3600

E: pmschina@pmeasuring.com

FRANCE

T: 33(0)1 60 10 32 96

E: pmsfrance@pmeasuring.com

GERMANY

T: +49 6151 6671 632

E: pmsgermany@pmeasuring.com

ITALY

T: +39 06 9053 0130

E: pmssrl@pmeasuring.com

JAPAN

T: +81 3 5298 8175

E: pmsjapan@pmeasuring.com

KOREA

T: +82 31 286 5790

E: pmskorea@pmeasuring.com

MEXICO

T: +52 55 2271 5106

E: pmsmexico@pmeasuring.com

All rights reserved.

SYSTEMS

PARTICLE

MEASURING

a Spectris company

T: +45 707 028 55

E: pmsnordic@pmeasuring.com

PUERTO RICO

T: +1 787 718 9096

E: pmspuertorico@pmeasuring.com

SINGAPORE

NORDIC

T: +65 6496 0330

E: pmssingapore@pmeasuring.com

SWITZERLAND

T: +41 71 987 01 01

E: pmsswitzerland@pmeasuring.com

TAIWAN

T: 886-3-5525300 Ext: 301

E: pmstaiwan@pmeasuring.com