

Microsoft Operating Systems KB and Patch Updates Release

Document ID#	PNP-6.0.0-Q4-2023
Software Name	PharmaceuticalNet PRO
Software Version	6.0.0
KB and patches release period	Q4 2023

Document Approval

Title	Name	Signature	Date
Product Line Manager	Daniele Pandolfi		February 27, 2024
Quality Representative	Anna Motti		February 27, 2024

Purpose

This document serves as a comprehensive release statement for Facility Pro Software, providing evidence of its compatibility with the latest Windows KB (Knowledge Base) and Patches. Compatibility is determined based on the assessment of the last three months' security patches listed in the official catalog of Microsoft. [Microsoft Update Catalog](#)

Through our testing, we have verified that Facility Pro Software seamlessly integrates with these updates, ensuring a smooth and uninterrupted user experience.

IMPORTANT NOTE

Official communication from Microsoft indicates possible challenges during installation of the Windows **KB5034441** update, causing overall update to fail or stuck.

For this reason, the **KB5034441** has been excluded from the Q4 2023 testing and shall not be installed.

List of tested KB and Patches

Server environment (Windows 10)

OS Name	Windows 10
OS Version	Version 22H2 (OS Build 19045.3930)
Hotfix(s)	Update KB5034466 Update KB5030841 Update KB4562830 Security Update KB4580325 Security Update KB4593175 Update KB5000736 Update KB5011048 Security Update KB5012170 Update KB5015684 Security Update KB5034122 Update KB5006753 Update KB5007273 Security Update KB5014032 Update KB5016705 Update KB5018506 Update KB5020372 Update KB5029709 Update KB5031539 Update KB5032907 Security Update KB5005699

Server environment (Windows Server 2016)

OS Name	Windows Server 2016 Datacenter
OS Version	Version 1607 (OS Build 14393.6614)
Hotfix(s)	Update KB5033910 Update KB3192137 Update KB4486129 Security Update KB4520724 Security Update KB4535680 Security Update KB4540723 Security Update KB4565912 Security Update KB4576750 Update KB4589210 Security Update KB5001402 Security Update KB5005698 Security Update KB5023788 Security Update KB5030504 Security Update KB5032391 Security Update KB5034119

Server environment (Windows Server 2019)

OS Name	Windows Server 2019
OS Version	Version 1809 (OS Build 17763.5329)
Hotfix(s)	Update KB5033911 Update KB4486153 Security Update KB4512577 Update KB4589208 Security Update KB5034127 Security Update KB5020374 Security Update KB5023789 Security Update KB5028316 Security Update KB5031589 Security Update KB5032306

RDP Server environment

OS Name	Windows Server 2019
OS Version	Version 1809 (OS Build 17763.5329)
Hotfix(s)	Update KB5033911 Update KB4486153 Security Update KB4512577 Update KB4589208 Security Update KB5034127 Security Update KB5020374 Security Update KB5023789 Security Update KB5028316 Security Update KB5031589 Security Update KB5032306

Client environment

OS Name	Windows 10
OS Version	Version 22H2 (OS Build 19045.3996)
Hotfix(s)	Update KB5034466 Update KB5028853 Update KB4562830 Security Update KB4580325 Security Update KB4593175 Update KB5011048 Security Update KB5012170 Update KB5015684 Update KB5034203 Update KB5014035 Update KB5023794 Update KB5028318 Update KB5031539 Update KB5034224 Security Update KB5005699

Conclusion

Particle Measuring Systems will maintain test evidence documentation on-site or at a secure off- site location and will require a non-disclosure agreement to be provided by those wanting access to the information.